


Camel Emergency Lifting Device



The Camel (Complete Air Moving Elevating Lift) is a state of the art air powered cushion designed to sit up and lift a fallen person from the floor.

Operating with low air pressure, the Camel offers an extremely stable and stressfree lift that can be performed either independently or with assistance from a carer.

With a maximum user weight of 320 kg, the Camel easily lifts people of any size and age but is particularly useful for lifting bariatric people, minimising the risks associated with moving and handling. Backrest provides extra support.

Specifications:

Product Code	2635
Max Height	560 mm
Min Height	Flat
Seat Depth	470 mm
Total Length	1120 mm
Width	700 mm
Weight	6.5 kg
SWL	320 kg