

Neos 150 Folding Hoist



Product Features:

The NEOS 150 Folding is a compact ergonomic designed medium sized heavy duty mobile patient lifting hoist with a safe working load of 150kg. This lifter is ideally suited as a compact and lightweight hoist and provides superior space. The swivel feature of the detachable yoke allows a patient to be rotated through 360° at full lift height, making the carer's job easier and safer.

Product Benefits:

- Quick release folding structure enables lift to be conveniently folded for storage and transfer
- Electric raise and lower powered by Linak
- Foot pedal control for leg adjustment
- 75mm front castors and 100mm locking rear castors
- Quick release yoke

Specifications:

Product Code	2878	Lift Range - Yoke	410 -1650 mm
Base Width Min	615 mm	Castor Size	75 / 100 mm
Base Width Max	965 mm	Under Bed Clearance	135 mm
Folded Overall Length	1230 mm	Lifter Weight	35 kg
Folded Overall Width	400 mm	SWL	150 kg