

Alto MKII 200 Hoist



**Heavy Duty
Patient Lifting Hoist.**

Alto MKII 200 Hoist

Product Features:

A powered pivoting frame makes this hoist the perfect choice for lifting larger patients. With the push of a button repositioning and tilting the patient is safe and strain free. A medium sized heavy duty mobile patient lifting hoist with a safe working load of 200kg and a compact ergonomic design. Alto™ Lifter whilst being compact provides superior patient lift space through the use of SPACELINK™ technology.


The Alto™ Lifter is the lifter of choice for institutional use or home care where a high capacity lift is required. Electric leg adjustment is a standard feature. It utilises the new SPACELINK™ concept, a breakthrough for optimised patient space, that enhances floor lift capability and increases max. lifting height. With maintenance free bushes in all the major moving parts this hoist is an Australian Made example of Quality and Strength that is hard to find in other hoists.



Standard Yoke



4 Point Yoke




Tri Hook Yoke



Pivot Frame



Powered Pivot Frame



Weigh Scale

Specifications:

Product Code	2891	Hoisting Range*	1620 mm
Base Length	1150 mm	Minimum Storage Height	1350 mm
External Base Width	645 mm	Under Base Height	50 mm
Maximum Internal Length	860 mm	Maximum Hoist Reach	1060 mm
Internal Base Width Closed	578 mm	Minimum Hoist Reach	610 mm
Internal Base Width Open	1130 mm	Max Hoist Reach	41 kg
Base Width Adjustment	Electric	Max Height	4400 mm
Under Bed Clearance	125 mm	Min Hoist Reach at Min Height	430 mm
Minimum Height o	420 mm	Handle Height	1000 mm
Spreader Bar	300 mm	Turning Radius	1210 mm
Maximum Height	870 mm	Mass of Heaviest Part (Base Assembly)	17.5 kg
Spreader Bar	1920 mm	Lifter Weight	42 kg
		SWL	200 kg