

# Forte 320 Lifting Hoist



**Heavy Duty  
Patient Lifting Hoist.**

# Forte 320 Lifting Hoist



Bariatric 4 Point Yoke



Featuring Linak Jumbo Care advanced diagnostic control system.



Bariatric 4 Point Yoke with Integrated Scale



SPACE LINK™  
Breakthrough for optimised patient space and enhances floor lift capacity

## Product Benefits:

With a safe working load of 320kg this extremely heavy duty mobile patient lifting hoist provides full bariatric capability whilst remaining a compact and ergonomic design. Manoeuvrability and exceptional under bed clearance are key design features.

FORTE 320 provides superior patient lift space through the combination of the recessed actuator and the use of SPACE LINK™ technology. The aesthetic design of the concealed electric leg mechanism provides infection control advantages.

FORTE 320 is the lifter of choice for institutional use where a high capacity bariatric lifter is required.

## Specifications:

<b>Product Code</b>	2898
<b>Base Length</b>	1420 mm
<b>External Base Width Closed</b>	715 mm
<b>External Base Width Open</b>	1200 mm
<b>Internal Base Width Closed</b>	600 mm
<b>Internal Base Width Open</b>	1080 mm
<b>Under Bed Clearance</b>	111 mm
<b>Minimum Yoke Height</b>	245 mm
<b>Maximum Yoke Height</b>	1940 mm
<b>Lifter Weight</b>	63 kg
<b>Minimum Storage Height</b>	1445 mm
<b>SWL</b>	320 kg