

# CuroCell SAM Air Mattress Overlay



Constructed from heavy duty aluminium, this bariatric seat walker can cater for those up to 220kg. It features a removable padded backrest and ergonomically shaped handgrips for user comfort. The handbrakes are lockable for stability. A wire basket is included also.

## PRODUCT BENEFITS:

- Suitable for high risk clients, and also medium risk as a preventative
- Dynamic air mattress overlay
- Treatment up to and including Pressure Injury Stage 2
- User weight 1kg to 200kg
- Comes complete with an anti-shear inner cover which has been proven to reduce pressure by up to 8%
- Easy to use dial pump with table to calculate most appropriate setting
- 2 year manufacturers warranty

## Specifications:

<b>Product Code</b>	3670
<b>Width</b>	850 mm
<b>Length</b>	2000 mm
<b>Thickness</b>	100 mm
<b>SWL</b>	200 kg

## Sizes Available:

<b>Product Code</b>	3671	3672	3673
<b>Width</b>	900 mm	1050 mm	1200 mm
<b>Length</b>	2000 mm	2000 mm	2000 mm