

# Regency Care Chair



**Patented Gel, Air &  
Foam Seating System**

# Regency Care Chair

The 5050G Regency Care Chair is designed to give the client comfort by combining the original tilt system with an independent multi-positional back rest. Its many features take into account the needs of both carer and user of the chair. The Regency 5050G derives its status from a proven and successful track record across a broad spectrum of the industry and has features such as:

- Forward transfer friendly
- Manoeuvrable caster base
- High quality upholstery
- Side Transfer Friendly
- Superior range of tilt and completely flat recliner



## Product Benefits:

- Reclining backrest – full recline to supine position
- Tilting seat
- Pressure management seating – gel cushion in seat, self-air adjusting back, profiled head cushion
- 2-position armrests
- Forward seat tilt for ease to stand, legrest angles under seat
- Elevating contoured legrest complete with fold out footplate
- Low skin shear 2-way stretch fabric on seating surfaces
- Drop down armrests and swingaway headwings – side transfer friendly
- Patient lifter friendly castor base
- 5 year limited warranty
- 175mm castors, braking on rear, directional lock on front
- Quality Australian made construction

## Specifications:

<b>Product Code</b>	5050G
<b>Seat Width</b>	460 mm
<b>Seat Depth</b>	520 mm
<b>Seat Height (Tilt/Normal)</b>	570 – 620 mm
<b>Backrest Height</b>	700 mm
<b>Armrest Height</b>	175 mm
<b>Legrest Length</b>	480 – 650 mm
<b>Overall Width</b>	660 mm
<b>Overall Depth</b>	1040 mm
<b>Overall Height</b>	1335 mm
<b>SWL</b>	180 kg

## Colours Available:



Mink



Lagoon



Mustard Pattern



Arcadia Marina



Green



Bouquet Navy



Purple



Tangerine