

Shopper Walker



The Shopper Walker is excellent for outdoor use, with its heavy duty frame and large 200 x 50mm puncture proof castors. Other features include large shopping basket, adjustable height handles, Orthopaedic hand grips, sturdy seat and cable brake lockable castors. A plastic tray is also included giving further versatility.

Product Benefits:

- Height adjustable ergonomic shaped handgrips
- Heavy duty powdercoated steel construction
- Lockable cable type brakes
- 200mm castors

Specifications:

Product Code	6740
Seat Height	600 mm
Overall Width	545 mm
SWL	100 kg