

# byACRE Overland Walker



**The sturdy, light  
weight and versatile  
multi-terrain walker.**

# byACRE Overland Walker

The byACRE Carbon Overland is a lightweight offroad walker to take you to the places you love. With the style of a Land Rover this walker is strong and durable and ready for your next adventure.

## Product Benefits:

- Lightest off-road walker at 6.7 kg for easy lifting in your vehicle
- Easy to fold, easy to grip and easy to lift
- Fully assembled, to make it hassle free to start exploring
- Extremely stable carbon fibre frame – 5 x stronger than steel
- Ultra-durable pneumatic rubber tires for better suspension
- Rigid front forks and sealed bearings
- Aircraft grade aluminium in gunmetal colour, for strength and lightness
- A safe working load of 120kg

## Specifications:

<b>Product Code</b>	6880LW
<b>Handle Height</b>	820 - 950 mm
<b>Length</b>	730 mm
<b>Seat Width</b>	470 mm
<b>Seat Height</b>	610 mm
<b>Width</b>	660 mm
<b>Width when Folded</b>	250 mm
<b>Wheel Diameter</b>	250 mm
<b>Weight</b>	6.7 kg
<b>SWL</b>	120 kg

## Also Available:

<b>Product Code</b>	6886BB	6886DG
<b>Colour</b>	Bentley Brown	Defender Green

