

Multicare X Bed



**A combination of
advanced innovations
& functional features.**

Multicare X Bed

A premium ICU bed, the Multicare X with unique frame-based lateral tilt. Accompanied by Opticare X, innovative constant low pressure supporting surface.

Product Benefits:

- Improve Pulmonary Outcomes: Frame-based lateral tilt with programmable Automatic Lateral Therapy, guided by Electrical Impedance Tomography supporting the reduction of complications associated with mechanical ventilation.
- Early Mobilization: Our Early Mobilization Program, featuring Mobi-Lift® handle or lateral tilt, can contribute to better and faster patient recovery.
- Pressure Injury Prevention: New generation of integrated mattress Opticare® X is fully automatic allowing the adjustment of mattress settings by the press of a single button to support individual patients clinical conditions and repositioning.
- Falls Prevention: Safety concept of the siderails, automatically engaged i-Brake® and multizone bed exit alarm keep the patient and caregivers safe in every situation.
- Improve Patient Experience: Features such as Ergoframe®, Opticare® X Mobile Mode and USB integration support patients physically and psychologically, improving their well-being and enhancing recovery.
- Lateral Tilt: The frame-based lateral tilt makes it easier and safer for nurses to provide daily care. It also reduces the physical strain applied on the caregiver's back, reducing risks of back injuries.
- X-ray translucent mattress platform
- Ergoframe seating system
- Multiprotect side rails with down sensor & tilt indicator
- Integrated weighing scales and bed exit alarm
- Multipurpose keyboard in both head siderails
- I-Power Drive Wheel
- Integrated USB ports



Specifications:

Product Code	74410
External Dimensions	2150 x 1050 mm
Maximum Safety Siderail Height (Above Mattress Platform)	450 mm
Bed Extension	220 mm
Mattress Dimensions	2140 x 875 mm
Maximum Mattress Height	230 mm
Underbed Clearance	150 mm
Height Adjustment	440 - 820 mm
Lateral Tilt	±30°
Thigh / Calf Rest Angle	30° / 38°
Backrest Angle	70°
TR / ATR Angle	13° / 16°
SWL	250 kg