

Bedside Commode



The standard bedside commode offers padded armrests for comfort and features adjustable height, a removable padded seat cover and a plastic toilet seat with bowl and lid.

Product Benefits:

- Stackable design
- Various colours available
- 125kg SWL

Specifications:

Product Code	7952PA
Overall Width	575 mm
Overall Height	940 - 1090 mm
Overall Depth	600 mm
Seat Width	440 mm
Seat Height	500 - 650 mm
Seat Depth	440 mm
Back Height	475 mm
Toilet Seat	475 - 625 mm
SWL	125 kg