

activeMEDICAL

Configura Comfort

Assessment Form



Form Information

Dealer Information

Date	/	/	Sales Rep	
Street Address			Email	
Suburb			Phone	
Postcode		State	Fax	

Client Information

Date	/	/	Full Name	
Contact			Street Address	
Funding Number			Suburb	
Plan Date - Start to Finish		to	Postcode	State
Plan Manager			Client Weight	kg

Funding Type

<input type="checkbox"/>	M.A.S.S.	<input type="checkbox"/>	NDIS	<input type="checkbox"/>	Homecare Package	<input type="checkbox"/>	Other: _____
--------------------------	----------	--------------------------	------	--------------------------	------------------	--------------------------	--------------

Prescriber Information

Date	/	/	Email	
Full Name			Street Address	
Organisation			Suburb	
Mobile			Postcode	State
Office Phone				

Comments

Specifications (Areas highlighted in bold are standard specification)

Chair Size	Seat Width	Part Number	Seat Depth	Seat Height
Small (130kg SWL)	457 mm (18")	5550S - Duratek	406 mm (16") ¹	388 mm (15") ⁴
			457 mm (18")	406 mm (16")
		5555S - Vinyl	457 mm (18")	431 mm (17") ²
			508 mm (20")	457 mm (18") ³
Medium (160kg SWL)	508 mm (20")	5550M - Duratek	406 mm (16") ¹	398 mm (17") ⁴
			457 mm (18")	457 mm (18")
		5555M - Vinyl	457 mm (18")	482 mm (19") ²
			508 mm (20")	508 mm (20") ³
Large (160kg SWL)	559 mm (22")	5550S - Duratek	406 mm (16") ¹	439 mm (17") ⁴
			457 mm (18")	457 mm (18")
		5555L - Vinyl	457 mm (18")	482 mm (19") ²
			508 mm (20")	508 mm (20") ³

Notes

Depth Adjustment ¹ - Kit 2" (50 mm) 5500 - 4

Height Adjustment ² - Kit 1" (25 mm) 5500 - 6

Height Adjustment ³ - Kit 2" (50 mm) 5500 - 7

Height Reduction ⁴ - Kit 1" (18 mm)

Pressure

Visco a Cushion (std)

Cushion Air + Pump CA2400

Integrated Self Supply



Posture

Waterfall Backrest (std)



Laterals CA2413



Come as each

Profiled Headrest



Lateral Support Backrest

Small
Medium
Large



Adjustable Lateral Support

Small
Medium
Large



Transfer Options

Description	Long Description	Part Number	Select
Drop-down Armside Left-hand side (as sitting in the chair)	Enables users to do a slideboard transfer in and out of the chair with ease. Can be bolted on at any time if the users needs change.	Small - Duratek Small - Vinyl Medium/Large - Duratek Medium/Large - Vinyl	<input type="checkbox"/> Small
			<input type="checkbox"/> Medium
			<input type="checkbox"/> Large
Drop-down Armside Right- hand side (as sitting in the chair)	Enables users to do a slideboard transfer in and out of the chair with ease. Can be bolted on at any time if the users needs change.	Small - Duratek Small - Vinyl Medium/Large - Duratek Medium/Large - Vinyl	<input type="checkbox"/> Small
			<input type="checkbox"/> Medium
			<input type="checkbox"/> Large

Extras

Armrests

Overhand on outside (std)

Overhand on inside - reduces width between
armrests by 1.5" (38 mm)

Rise Setting

Rise & Tilt **(std)** 23° (160kg)

Medium Rise 14° (135kg*)

Flat Rise 2° (110kg*)
* Medium & Large Chairs

Cognitive

Legrest Setting

Standard

Standard + 1" (25 mm)

Standard + 2" (50 mm)

2 - Button Handset



Support

Legrest Channel

Small

Medium

Large



Handset Control

Right Hand Side **(Std)**

Left Hand Side

Carer Position

Chair & Delivery Environmental Assessment

When carrying out an assessment please fill out the below form prior to ordering the chair.

Please give as much detail as possible.

Access

Is there parking / access for a large van?

Are there steps to the property?

Is it easy access to the property?

Internal

Are the doors and corridors wide enough for a chair to go through? (Min. width for most chairs is 650mm.)

Does the chair have to go upstairs? If yes, please specify number of steps and if a stair lift is fitted.

Is space being cleared for the chair?

Is there a suitable power supply for the chair?

Yes

No

Infection Control

Are there pets?

Yes

No

Is the house smoke free?

Yes

No

Other Comments

Save As

Print

Email

activeMEDICAL
SUPPLIES > SOLUTIONS > SUPPORT

P. 1800 267 267 | **F.** 07 3899 8825
E. hello@activemedicalsupplies.com.au

www.activemedicalsupplies.com.au