

# Action X-act Lite Cushion



Action Products XACT Cushion range continues to expand with the addition of the XACT LITE, a cushion for the active and mobile patient, including those who require sports-friendly cushions, and can self-propel, transfer and shift. The new design of the AKTON polymer reduces shear and pressure and provides maximum surface flexibility with no hammocking effect.

Additional features include light leg troughing which improves stability and posture without interfering with transfers, low profile design which helps foot propulsion, pre-ischial ridge to help maintain pelvic neutrality, dual-density foam base with firm rear corners to improve stability, and soft front corners to aid transfers. Cushion weight of 1.45kg, and height of 5.4 cm.

## PRODUCT BENEFITS:

- 54mm thick
- Suitable for clients medium - high risk
- 114kg weight capacity
- Easy to care for
- BodySmart cover included

## AVAILABLE SIZES:

16" x 16"	3580
18" x 16"	3581
18" x 18"	3582
20" x 18"	3583

**activeMEDICAL**

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au  
activemedicalsupplies.com.au