


Allegro Pivot Frame



The pivot frame is used with the pivot sling. The frame spreads the sling attachment points to ensure the patient's hips are not flexed excessively, avoiding the discomfort some patients experience with standard slings. The ability to position a person in a sitting or lying position, greatly assists in providing correct seating and positioning.

activeMEDICAL

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au
activemedicalsupplies.com.au