

Trust Care Outdoor Walker



It is light and elegant, both folded and unfolded. It has an X frame cross folding system which makes it easy to fold, takes up little space and stands alone when folded. Big front and rear wheels make it easy to get over curbs and other obstacles. Even more important is the spring system in the wheels that gives a smooth movement compensating for stiff joints. The walker has soft, comfortable rubber handles that are adjustable in height from 78 to 96 cms to fit different heights.

It has an easy to manage brake which of course can be locked and a comfortable seat in imitation leather for strength. This walker has an elegant design in every aspect.

Available in Black/Silver (60120) or Mother-of-Pearl/Silver – pictured (60150).

PRODUCT SPECIFICATIONS:

TECHNICAL DATA	
Product Code	60150
Seat Width	400mm
Seat Height	600mm
Handle Height	780mm – 960mm
Walker Weight	6.5kg
SWL	100kg

activeMEDICAL

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au
activemedicalsupplies.com.au