

# Sertain S4400 Series Bariatric Care Chair



For the bariatric client, this model will support up to 300kg, with functionality, comfort and the needs of both user and carer in mind. The electric operation of hi-lo adjustments mean transfers and patient care are able to be carried out at the correct working height by staff.





### FEATURES AND BENEFITS:

- Electric hi-lo 540-940mm – low enough to stand from with ease, but high enough to achieve a comfortable working height for nursing staff
- Electric reclining backrest – full recline for supine position
- Electric tilting seat
- Electric elevating legrest
- Handset operated with preset auto chair and auto layflat functions
- Pressure management seating – gel cushion in seat, air in backrest
- Fold out, extendible footplate
- Low skin shear 2-way stretch fabric on seating surfaces
- Removable two-position armrests
- Swingaway three position headwings
- Adjustable head rest cushion
- 150mm castors, 3 x braking, 1 x directional lock
- Battery back up system
- Quality Australian made product
- Warranty – 5 year Structural Frame, 2 years Castors, 3 years Electric System

### PRODUCT SPECIFICATIONS:

TECHNICAL DATA	
Product Code	5135
Seat Width	580mm
Seat Depth	520mm
Seat Height (Tilt/Normal)	550mm / 950mm
Backrest Height	700mm
Armrest Height	160mm
Legrest Length	440mm
Overall Width	790mm
Overall Depth	1080mm
Overall Height	1335mm
SWL	300kg

### OPTIONS:

- Layflat Support Bar
- Droparm
- Flip Up Tray Arm
- Oxygen Bottle Carrier
- IV Pole and Bracket



Layflat Support Bar



Flip Up Tray Arm



Droparm Option