

Silicone Fibre Elbow / Heel Protector



3105



3155

FEATURES & BENEFITS:

- Heel protector also available - 3155
- Provides comfort and assists in the prevention of pressure build up in the elbow area
- Adjustable Velcro securing straps to suit
- One size fits all Machine washable

This elbow cushion is made of millions of silicone coated cylindrical hollow cored fibres. These resilient fibres provide excellent pressure relief while the silicone coating acts as a lubricant to eliminate friction and reduce shear forces. The internal seams are carefully positioned to avoid creating pressure points.

activeMEDICAL

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au
activemedicalsupplies.com.au