

Fortress 5 Section Table



5 section support surface with electrical height adjustment for unrivalled versatility

- Drop down arm and shoulder rests allows a more comfortable position for patient and enables the therapist to get closer to client and reduce back strain
- Fitted with 4 retractable castors and central lifting system
- Head section with ergonomically designed removable cushion and facial opening
- Nicely curved cushions with premium padding
- Durable, anti-bacterial, washable and flame retardant upholstery
- Excellent stability with new 'Safetyline' streamlined base.

PRODUCT SPECIFICATIONS:

TECHNICAL DATA	
Table Width	630mm
Head Tilt Angle	-50° - 87°
Min Height	495mm
Max Height	940mm
Head Section Length	650mm
Mid Section Length	680mm
Leg Section Length	660mm
Overall Width	700mm
Overall Length	2060mm
SWL	225kg

activeMEDICAL

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au
activemedicalsupplies.com.au