

Regency Care Chair R3000 Series



PRODUCT BENEFITS:

- Reclining backrest – full recline
- Tilting seat
- Generous 490mm wide seat width
- Pressure management seating – gel cushion in seat, self-air adjusting back, profiled head cushion
- 2-position armrests
- Forward seat tilt for ease to stand, legrest angles under seat
- Elevating contoured legrest complete with fold out, extendible footplate
- Low skin shear 2-way stretch fabric on seating surfaces
- Drop down armrests and swingaway headwings – side transfer friendly
- Covers replaceable
- Patient lifter friendly castor base
- 5 year limited warranty
- 175mm front castors that are directional lock
- Quality Australian made construction

The R3000 Series Regency Care Chair is designed to be carer and client friendly by incorporating the original tilt flotation seat system and infinite reclining backrest into a multi-positional, multi-functional and manoeuvrable care chair. Large rear wheels make this model ideal for outdoors.

PRODUCT SPECIFICATIONS:

TECHNICAL DATA	
Product Code	5060
Seat Width	490mm
Seat Depth	520mm
Seat Height (Tilt/Normal)	570mm / 630mm
Backrest Height	700mm
Armrest Height	175mm
Legrest Length	480mm-650mm
Overall Width	660mm
Overall Depth	1000mm
Overall Height	1335mm
SWL	180kg

activeMEDICAL

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au
activemedicalsupplies.com.au