

Sliding Transfer Bench



PRODUCT BENEFITS:

- Plastic moulded textured top features holes to drain water away from the seating area.
- Locking device to stop seat movement at either end
- Velcro lap belt which provides increased security
- Curved backrest provides added support
- Integrated grab rail (specify left or right side)

The Sliding Transfer Bench provides the solution for individuals with restricted ability to use a standard transfer bench. The sliding seat allows easy transfer to the bathing position and back again. It has a height adjustable aluminium frame to suit individual user requirements.

SPECIFICATIONS:

TECHINICAL DATA

Product Code	1016
Seat Width	480mm
Seat Depth	330mm
Seat Height	485mm - 600mm
Bath Rim Clear	390mm - 505mm
Seat Back Height	380mm
Overall Width	1000mm
Overall Depth	390mm
SWL	150kg

activeMEDICAL

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au
activemedicalsupplies.com.au