

Commode Transfer Bench



PRODUCT BENEFITS:

- Two piece padded seat allows water to drain away from seating area
- Padded seat for user comfort
- Front opening for easier personal cleaning
- Rack that allows a bowl to slide under for commode use

The Commode Transfer Bench has all the features of the standard transfer bench with the extra benefit of padded seat with open front.

SPECIFICATIONS:

TECHINICAL DATA	
Product Code	1019
Seat Width	630mm
Seat Depth	405mm
Seat Height	460mm - 560mm
Bath Rim Clear	340mm - 440mm
Seat-Back Height	380mm
Overall Width	730mm
Overall Depth	600mm
SWL	125kg

activeMEDICAL

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au
activemedicalsupplies.com.au