

Aluminium Shower Chair with Arms



The ACare Shower Chair provides a comfortable seat for those unable to stand whilst showering. The armrests provide security when entering or exiting the shower, and assist the user on and off the seat.

- Seat height and angle adjusts to suit the individuals requirements
- Aluminium construction for long life
- One piece moulded supportive seat and backrest
- Small drainage slots allow water to escape and remove risk of entrapment
- Chair design is stackable for space saving

Specifications:

TECHNICAL DATA

Product Code	1315
Seat Width	440 mm
Seat Depth	430 mm
Seat Height	400 - 550 mm
Width Between Arms	460 mm
Backrest Height	350 mm
Overall Width	570 mm
Overall Height	540 mm
SWL	200 kg

> This product has been independently tested according to the Australian and New Zealand Standards for Shower/Toilet Chairs and Stools (Static): AS/NZS 3973:2009, AS/NZS 4688.2-2000, AS/NZS 4688.3-2000 and carries a five year warranty on the effectiveness of preventing the development of rust.