

# Extra Care Shower Chair with Arms



## PRODUCT SPECIFICATIONS:

TECHNICAL DATA	
Seat Width	440mm
Seat Depth	430mm
Seat Height	400mm - 550mm
Width Between Arms	500mm
Backrest Height	350mm
Overall Width	620mm
Overall Depth	570mm
SWL	135kg

The ACare Shower Chair provides a comfortable seat for those unable to stand whilst showering. The armrests provide security when entering or exiting the shower, and assist the user on and off the seat.

- Strengthened frame for users needing extra room and weight capacity
- Seat height and angle adjusts to suit the individuals requirements
- One piece moulded supportive seat and backrest
- Small drainage slots allow water to escape and remove risk of entrapment

**activeMEDICAL**

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au  
activemedicalsupplies.com.au