

Aluminium Shower Stool



The ACare shower stool with arms, provides security when entering or exiting the shower and assist the user on and off the seat.

- Compact size for the smaller shower recess
- Padded seat for user comfort
- Adjustable seat height for the individual's requirements
- Aluminium frame for longevity

Also available:

1551 PLASTIC SEAT OPTION

PRODUCT SPECIFICATIONS:

TECHNICAL DATA	
Product Code	1567
Seat Width	370mm
Seat Depth	290mm
Seat Height	520 - 670mm
Width Between Arms	460mm
Armrest Height (from floor)	640-840mm
Overall Width	550mm
Overall Depth	370mm
SWL	125kg

activeMEDICAL

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au
activemedicalsupplies.com.au