

# Shower Stool with Arms and Back



The ACare shower stool with arms and backrest, provides security when entering or exiting the shower and assist the user on and off the seat.

- Compact size for the smaller shower recess
- Padded seat and backrest for user comfort
- Adjustable seat height for the individual's requirements

## PRODUCT SPECIFICATIONS:

TECHNICAL DATA	
Product Code	1569
Seat Width	370mm
Seat Depth	290mm
Seat Height	520mm - 670mm
Width Between Arms	640mm - 840mm
Armrest Height	640mm - 840mm
Overall Width	570mm
Overall Depth	450mm
SWL	125kg

**activeMEDICAL**

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au  
activemedicalsupplies.com.au