

# Aluminium Over Toilet Frame



This aluminium over toilet frame is adjustable in height and fits over the toilet. This frame has been designed to assist the user who has difficulty sitting down or standing up by raising the toilet seat height.

## Product Benefits:

- Aluminium construction for long life
- No risk of rust
- Removable splash guard included
- Lift up seat that removes for cleaning
- Ideal for those with hip or knee injuries or who require an increased seat height
- Unique stacking design saves space when storing more than one frame
- Multi D leg adjustment ensures stability on uneven surfaces
- Push button height adjustment
- Non slip guides on base of legs

## Specifications:

### TECHNICAL DATA

Product Code	1755
Width Between Arms	450 mm
Seat Height	475 - 625 mm
Overall Width	560 mm
Armrest Height	700 - 850 mm
Footprint (WxD)	560 x 590 mm
SWL	150 kg