

# Juvo Transfer Bench



The Juvo Heavy Duty Transfer bench is both compact and sturdy providing safe and easy access for bathing or showering.

With a weight limit of 200kg this suits a vast range of clients whilst remaining lightweight for easy manoeuvrability. The Heavy Duty bench features height adjustable legs, non slip feet, the option of suction cups and reversible backrest to suit all bath types. With a textured non-slip seat with drainage holes and side rail the transfer bench provides stable seating for our clients.

## SPECIFICATIONS:

TECHINICAL DATA	
Product Code	1785
Seat W x D	680 x 410mm
Seat Height	470 - 570mm
Backrest Height from Seat	360mm
Armrest Height from Seat	220mm
Bath Rim Clearance	390-490mm
Overall W x D	780 x 550mm
Overall Height	930mm
SWL	200kg

**activeMEDICAL**

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

mfs@activemedicalsupplies.com.au  
activemedicalsupplies.com.au