

# Throne Arms



The Throne is a 3 in 1 toilet support rail system, providing assistance to stand from the toilet. The rails can be positioned in the standard, splayed or fold down position, making this a versatile and adaptable piece of equipment. It is easily attached to the toilet utilising the existing toilet seat fixings.

## FEATURES & BENEFITS:

- Durable powder coated steel construction
- Ergonomic arm design for maximum support
- Rubber lined clamps for stable fitting
- Throne spacer available in 50 or 80mm height if increased toilet seat height is necessary.

## PRODUCT SPECIFICATIONS:

TECHNICAL DATA	
Product Code	1830
Arm Height	260mm
Width Between Arms (Upright / Splayed)	510mm / 640mm
SWL	110kg

**activeMEDICAL**

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au  
activemedicalsupplies.com.au