



ActiveCare Pivot Frame Sling

ActiveCare Mesh

Extra Extra Small	227239
Extra Small	227240
Small	227241
Medium	227242
Large	227243
Extra Large	227244
Extra Extra Large	227238

ActiveCare Fabric

Extra Small	227130
Small	227133
Medium	227134
Large	227135
Extra Large	227136
Extra Extra Large	227132



The Pivot sling is designed for use with the pivot/tilt frame. It attaches with the keyhole tags engineered to the highest standards in strength and size tolerances and feature auto clip on action (no safety straps required). It also features sling centre marking and positioning handles to assist correct positioning and woven label that will not wash out or fade. Available in sizes from x-small through to XX-large and in fabric or mesh styles.

	XXS	XS	S	M	L	XL	XXL
A	697	797	897	997	1120	1240	1360
B	170	190	210	230	245	250	250
C	250	685	790	895	995	1080	1125
D	875	1060	1230	1385	1525	1650	1725
E	295	375	440	490	530	570	600

activeMEDICAL

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au

activemedicalsupplies.com.au