


ActiveCare Hygiene General Purpose Head Support Sling




ActiveCare Mesh

Extra Small	227405
Small	227406
Medium	227358
Large	227359
Extra Large	227407
Extra Extra Large	227408

ActiveCare Fabric

Extra Small	227401
Small	227402
Medium	227395
Large	227396
Extra Large	227403
Extra Extra Large	227404



The ActiveCare General Purpose Hygiene sling has been designed for use with patients who need to have the sling remain in place whilst being toileted. The sling has head support and is designed so that clothing can be easily removed for toileting. The Hygiene model comes standard with a torso strap for optimum safety. The loop style attachments mean this is a sling used on a patient lifter with a yoke style attachment. The multiple loops allow for correct positioning of the client.

	XS	S	M	L	XL	XXL
A	800	900	1000	1120	1240	1360
B	195	205	215	225	235	245
C	500	585	670	740	810	850
D	1020	1185	1350	1520	1650	1720
E	520	600	680	780	840	870