



ActiveCare Universal Deluxe Sling with Head Support



This sling is good for all standard transfers and is exceptionally comfortable for the patient. The leg sections are shaped to provide good support to the hip and thigh areas. It is designed to hold the leg more securely than standard sling leg sections, as the leg is held from both sides instead of just the inside of the leg. The leg sections also spread the load over a larger area on the patient's legs so tend to reduce the pressure on the underside of the legs.

	XS	S	M	L	XL	XXL
A	715	800	900	1000	1105	1205
B	160	185	210	235	260	285
C	685	815	945	1075	1205	1270
D	1000	1175	1355	1535	1710	1820
E	295	340	385	430	475	520
F	375	440	500	560	625	690



ActiveCare Mesh

X-Small	227430
Small	227431
Medium	227432
Large	227433
X-Large	227434

ActiveCare Fabric

X-Small	227435
Small	227436
Medium	227437
Large	227438
X-Large	227439

ActiveCare Parasilk

X-Small	227440
Small	227441
Medium	227442
Large	227443
X-Large	227444

activeMEDICAL

SUPPLIES > SOLUTIONS > SUPPORT