


Allegro Salsa Sling



The Salsa sling is for use with the Salsa stand up hoist. It features optimum back support and under arm padding for a safe and comfortable lift and also allows for the easy removal of clothing for toilet transfers.


The Allegro Salsa sling features a torso band for security and stability with buckle (b) or Velcro (v) fastening

Allegro Buckle

Small	27701
Medium	27702
Large	27703
Extra Large	27704

Allegro Velcro

Medium	27705
Large	27706
Extra Large	27707



	M	L	XL
A	880	1060	1220
B	290	370	410
C (b)	960	1020	1220
C (v)	1180	1320	1440

activeMEDICAL

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au
activemedicalsupplies.com.au