

Sonata Hoist



A lightweight, compact mobile patient lifting hoist with innovative and ergonomic design the Sonata 150 is an easy to use everyday lifter. With a safe working load of 150 kg the Sonata offers an outstanding lift capacity in a small patient lifter. It utilises the SPACELINK™ concept, a breakthrough for optimised patient space, that enhances floor lift capability and increases maximum lifting height.

Superb value. Even better quality.

With maintenance free bushes in all the major moving parts this hoist is an Australian Made example of quality and strength that is hard to find in other patient lifters.

The Allegro range of patient lifting equipment including lifters and accessories has been designed and tested to comply with:
AS ISO 10535-2002 Hoists for the transfer of disabled persons – Requirements and test methods.

- Manual leg adjustment is a standard feature.
- Supplied complete with battery, in-built battery charger and instruction manual.
- 150kg SWL

Email us for an order or quote enquiries@activemedicalsupplies.com.au

Compatible with powered pivot, standard yoke and manual pivot frame (pivot and yoke not included)



Low profile castor assist hoist to fit under low beds that are required by many patients



Interactive controls that show how much battery is left.



SpaceLink increases maximum lifting/height



PRODUCT SPECIFICATIONS:

TECHNICAL DATA	
Product Code	2880
Base Length	1005mm
External Base Width Closed	605mm
External Base Width Open	1070mm
Internal Base Width Closed	510mm
Internal Base Width Open	975mm
Base Width Adjustment	3 Position Manual
Under Bed Clearance	130mm
Minimum Height of Yoke	420mm
Maximum Height of Yoke	1735mm
Hoisting Range*	870mm
Minimum Storage Height	1260mm
Under Base Height	55mm
Maximum Hoist Reach	875mm
Minimum Hoist Reach	555mm
Hoist Reach at Maximum Height	555mm
Mass of Heaviest Part (Mast Assembly)	22kg
Lifter Weight	41kg
SWL	150kg

activeMEDICAL

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au
activemedicalsupplies.com.au