

# FORTE 320 Hoist



Featuring Linak Jumbo Care advanced diagnostic control system.

## Heavy Duty Patient Lifting Hoist

# Features and Specifications



## BENEFITS & USAGE:

With a safe working load of 320kg this extremely heavy duty mobile patient lifting hoist provides full bariatric capability whilst remaining a compact and ergonomic design. Manoeuvrability and exceptional under bed clearance are key design features.

FORTE 320 provides superior patient lift space through the combination of the recessed actuator and the use of SPACE LINK™ technology. The aesthetic design of the concealed electric leg mechanism provides infection control advantages.

FORTE 320 is the lifter of choice for institutional use where a high capacity bariatric lifter is required.



Bariatric 4 Point Yoke




Bariatric 4 Point Yoke with Integrated Scale



2 Point Yoke with Integrated Scale



Concealed Electric Leg Spread Mechanism



Low Profile Leg – Twin 75mm Castors



### SPACE LINK™

- breakthrough for optimized patient space
- enhances floor lift capacity

## PRODUCT SPECIFICATIONS:

### TECHNICAL DATA

	100mm Front Castor	Twin 75mm Castors
Base Length	1433mm	1420mm
External Base Width Closed	715mm	715mm
External Base Width Open	1205mm	1200mm
Internal Base Width Closed	590mm	600mm
Internal Base Width Open	1080mm	1080mm
Under Bed Clearance	126mm	111mm
Minimum Yoke Height	260mm	245mm
Maximum Yoke Height	1955mm	1940mm
Lifter Weight	63.5kg	63.0kg
Minimum Storage Height	1460mm	1445mm
SWL	320kg	320kg

**activeMEDICAL**

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au  
activemedicalsupplies.com.au