

staebel
Chair Scale



The Staebel Chair Scale with its patented design maintains reliability and uncompromising precision. All of the user's care needs, including safety, comfort, caregiver's use and maintenance, are met. The compact design ensures the most cost efficient shipping method and reduces the risk of carer strain particularly when weighing an obese/bariatric patient.

Unique base frame design gives care staff assurance that scale is compatible with other equipment utilised in care setting (e.g.. Full-sling/Standing hoist, Regency chair)

Product Benefits:

- Twin wheel castors are standard and do not create excessive noise when wheeling scale along a corridor
- One piece footplate flips up and has been designed with skin tear prone patients in mind
- Armrests flip back to allow scale to be utilise in bariatric care setting
- Base frame comes standard with a level indicator
- BMI function/ Tare/Zero/Hold/Print/kg/lb
- Scroll Down/Up on LCD screen
- 5 digit ¾" LCD screen that can be tilted to accommodate different care staff heights
- Automatic re-weigh
- Up-datable firmware
- Ports: RS232 and USB
- Amazing lightweight yet durable construction - only weighs 17.5 kg
- 24 month warranty
- 300kg SWL