

Action Gel Pilot Cushion



The Pilot is a low-profile cushion that provides exceptional pressure relief for most users with a pressure ulcer history. 2.54 cm of Akton polymer provides immersion for bony prominences and protects clients with fragile skin integrity. This multi-use cushion is recommended for use on manual and power wheelchairs, truck or taxi seats, or for other uses requiring significantly reduced high-pressure areas, vibration, and shear. The low profile allows users to effectively reach the floor for more efficient foot-assisted propulsion. It also conforms to pre-contoured surfaces to maximize stability and positioning when required.

PRODUCT BENEFITS:

- 1" (25mm) thick
- Suitable for medium-high risk clients
- Not restricted by any weight limit
- Easy to care for
- Quality, durable Pilot cover included. Cover can be machine washed (incontinence covers are an optional extra)

SIZES AVAILABLE:

- 16x16" 3530
- 18x16" 3531
- 18x18" 3532
- 20x18" 3533
- Custom sizes also available.

activeMEDICAL

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au
activemedicalsupplies.com.au