


# Relax Trio Cushion





The Relax Trio is an ultralight cushion that benefits from four sections of FloGel being positioned into a polyurethane base, with a viscoelastic foam upper layer. It is specially designed to conform to the shape of the user, increasing surface area and relieving pressure. The cushion can be rotated, it has no designated front or back.


## PRODUCT BENEFITS:

- 1.5" (38mm) thick
- Suitable for medium clients
- Clear protective PU film coating is easy to wipe clean and is suitable for incontinent users
- Lightweight and easy to transport
- Includes outer cover
- 80kg SWL

# Features and Specifications



Pressure diagram with Relax Trio cushion



Pressure diagram without wheelchair cushion



## RANGE OF SIZES AVAILABLE:

Product Code	Dimensions	Cushion Weight
35330	40cm x 40cm x 4cm / 16" x 16" x 1.5"	1.3kg
35332	45cm x 40cm x 4cm / 18" x 16" x 1.5"	1.5kg
35333	45cm x 45cm x 4cm / 18" x 18" x 1.5"	1.6kg