

# Relax Duo Profile Cushion Cushion



The Duo-Profile offers a pre-contoured surface giving excellent hip and leg alignment.

Has a 3 Layered Design:




- The base is Fast Recovery Foam
- The mid-section is Slow Foam/ Memory Foam
- The top layer is Medical Grade Silicone

The Duo-Profile uses a collapsible gel topper. This unique design allows you to quickly immerse into the cushion, while it conforms to your shape.

## PRODUCT BENEFITS:

- 3 Layered Design
- Includes outer cover
- 130kg SWL

# Features and Specifications



Integrated collapsible gel structure.

The combination of foam & gel can be effective in preventing both skin changes and pressure sore formation.



## RANGE OF SIZES AVAILABLE:

Product Code	Dimensions	Cushion Weight
35380	40cm x 40cm x 8cm / 16" x 16" x 3"	1.8kg
35382	45cm x 40cm x 8cm / 18" x 16" x 3"	1.9kg
35383	45cm x 45cm x 8cm / 18" x 18" x 3"	2.4kg