


Action Gel Commuter Cushion



PRODUCT BENEFITS:

- 57mm thick
- Suitable for medium risk clients
- Not restricted by any weight limit
- Easy to care for
- Basic cover included

The Action Commuter offers basic pressure and shear reduction for low risk users. Combined 1.27cm Akton polymer and 4.45cm foam are sealed in a water proof film for incontinence protection.

OTHER SIZES AVAILABLE:

16" x 16"	3560
18" x 16"	3561
18" x 18"	3562
20" x 18"	3563
20" x 20"	3564

activeMEDICAL

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au
activemedicalsupplies.com.au