

SAM BAS PRO Mattress & Pump



Mattress Replacement System
with Very High Comfort.



Features and Specifications

BENEFITS & USAGE:

The CuroCell S.A.M.® PRO mattress replacement system with an integrated foam core can be used as an effective aid in the prevention and in the treatment of pressure ulcers up to and including category 3. It can also be used as an aid for pain therapy.

This constant low pressure mattress offers an exceptionally high level of comfort as well as excellent pressure redistribution. With the vastly improved higher maximum weight capacity of 250 kg, this overlay will meet the ever increasing challenges of heavier patients.

The CuroCell S.A.M.® PRO features a new cost effective pump that includes both CPR and Transport functions. The pump automatically determines when to change pressure in contrast to running constantly.

Features and Characteristics

- CPR function – rapid deflation
- Soft, elastic, vapour permeable hygienic cover
- Very good cleaning properties
- Silent running and durable pump
- Easy to handle and maintain
- Unique labelling
- Energy efficient
- Maximum Pressure (caring mode)
- Extra soft inner cover
- Heel function

Facts:

- Type: Constant low pressure mattress
- Prevention and treatment up to and including category 3
- Warranty: 2 years
- Cleaning of cover: Wipe with cleaning agent and/or disinfectants. Machine wash max 95 °C, tumble drying.

PRODUCT SPECIFICATIONS:

TECHNICAL DATA	
Product Code	3683
Height	160mm
Power Supply	230V / 50Hz
SWL	250kg

SIZES AVAILABLE:

800 x 2000mm	800 x 2100mm
850 x 2000mm	850 x 2100mm
900 x 2000mm	900 x 2100mm
1050 x 2000mm	
1200 x 2000mm	

OPTIONAL ACCESSORIES

- Durable bag for transportation and storage
- Hygiene cover with welded seams

activeMEDICAL

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au
activemedicalsupplies.com.au