


# Positioning and Support Wedge



## FEATURES AND BENEFITS:

This foam positioning support wedge is used to place a user in a side position in order to relieve the sacrum, which is an extra sensitive area. For the side position, two support wedges are used together. The support wedge is easy to use and easy to clean. It is supplied with a zipped, liquid resistant and vapour permeable cover. The cover can be wiped or machine washed.

- 70cm x 26cm x 15-5cm
- Soft comfortable foam, stretch vinyl covered



**activeMEDICAL**

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au  
activemedicalsupplies.com.au