


Visko Mattress



The Viskomatt includes a double layer foam mattress including a 40 mm top layer of visco elastic foam, suitable for those spending a long time in bed, providing optimum pressure distribution across all body area in contact with the mattress.

A wide version is also available (3755W) with overall dimensions of 2000mm long x 1000mm wide x 140mm high.

PRODUCT BENEFITS:

- Suitable for medium risk patients
- Two way stretch PU waterproof cover
- Full length flap conceals the zip for infection control reasons
- Fully removable zippered cover
- Mattress cover is washable (up to 95 degrees C) and can be tumble dried. Moisture protection properties mean the cover has a high resistance to bacteria, viruses and fungi
- Mattress does not require turning
- 12 month warranty
- Overall size: 1990mm long x 870mm wide x 140mm high

activeMEDICAL

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au
activemedicalsupplies.com.au