

# Relax Tilt and Recline



The Active Care Relax is the latest in modular design. This quality aluminium wheelchair has a gas assisted tilt and recline adjustment built into the one product. The many adjustments within the chair features means it can be adapted to the individuals needs, removing the need to opt for a custom made product in many cases.

- Gas assisted seat tilt and recline
- Contoured adjustable headrest
- Height adjustable comfort backrest with lateral support
- Depth adjustable comfort seat with pommel
- Multiple axle positions provide adjustable seat height and wheelbase length
- Quick release rear wheels
- Adjustable height armrests with padded side guards
- wingaway removable elevating legrests with calf pads and angle adjustable footplates
- Adjustable height push handles
- Anti tip wheels
- 140kg capacity
- \*Also available in 51cm seat width

**activeMEDICAL**

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au  
activemedicalsupplies.com.au