

Embark Bariatric Sitting Chair



Heavy Duty
Patient Sitting Chair

Features and Specifications



Contoured Backrest angle adjustment



Contoured seat



Mobility wheels on back legs



Drop down, removable arms for side transfer, adjustable seat



Rear push bar



Adjustable height legs

PRODUCT SPECIFICATIONS:

TECHNICAL DATA	
Product Code	5966
Seat Width	780mm
Seat Height	460-540mm
Seat Depth	600mm
Back Height	675mm
Overall Width	900mm
Overall Height	100-1085mm
Width Between Arms	790mm
SWL	425kg

BENEFITS & USAGE:

Specially developed for hospital use, with a strong frame and high weight capacity, vinyl upholstery and tilt and move wheels to cater for the bariatric patient.

- Adjustable height
- High Contoured backrest incorporating lumbar support
- Contoured seat
- Padded drop down, removable arms for side transfer
- Adjustable seat depth
- Standard with tilt and move wheels for ease of manoeuvrability