

# Unilite Forearm Walker



## FEATURES & BENEFITS:

- Height adjustable ergonomic shaped hand grips
- Powder coated steel construction
- Lockable cable type brakes
- 200mm castors

## PRODUCT SPECIFICATIONS:

### TECHNICAL DATA

Product Code	6742
Seat Width	400mm
Seat Height	600mm
Handle Height	890-1100mm
Walker Weight	14kg
SWL	100kg

The ACare Unilite Forearm Walker is designed to support and tutor the user in an upright gait position. The very stable base gives added confidence and security to the user. Features include padded gutter supports that are adjustable in the forearm length, height and width. The popular folding design allows the Unilite to fold easily for transport and storage. This model is fitted with four large wheels and lockable cable brakes.

**activeMEDICAL**

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au  
activemedicalsupplies.com.au