

# Navigator Forearm Walker



## SPECIFICATIONS

Seat Width	450mm
Overall Width	550mm
Overall Depth	660mm
Seat Height	620mm
Handle Height	950-1150mm
Walker Weight	10.8kg
Overall Width when folded	310mm
SWL	150kg

- The navigator is foldable, staying firmly closed making it is easy to transport.
- The soft and comfortable arm supports and the brake handles are multi adjustable in depth, angle and height. The frame and seat are tested and approved on users with a maximum weight of 150 kg.
- Light and reliable brakes act on the TPE covered unbreakable wheels.
- The optimal and rigid design of the frame and the fore arm support ensures the best possible user comfort (lightness, turning circle and steadiness).
- The Navigator is easily height adjustable in the range from 90 – 115 cm.

**activeMEDICAL**

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au  
activemedicalsupplies.com.au