

Server HD Rollator



SPECIFICATIONS

Overall Measurements: (width x depth x height)	69 x 68 x 74–102 cm
Seat width	54 cm
Transportation measurements folded	Depth 68cm Width 24 cm
Weight	From 7.9 kg
Shopping bag: (width x depth x height)	35 x 17 x 25 cm
Max load	5kg
Max. user weight	200 kg
Colour	Sober Grey

Server HD folds easily and stays that way with a simple locking system that doubles up as an ergonomically shaped carrying handle.

- Frame and seat are tested and approved for maximum user weight of 200 kg.
- The brake mechanism is light yet positive. TPE-covered, unbreakable wheels.
- Server HD push handles are height adjustable with 10 positions, from 74 cm to 102 cm.
- The seat height is 62 cm with an extra wide 54 cm seat base. User comfort is ensured with an optimal frame design and ergonomic push handles.

activeMEDICAL

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au
activemedicalsupplies.com.au