

# Ligerio Seat Walker 55cm



The Ligerio a popular seat walker manufactured from durable lightweight aluminium. The larger castors make this an ideal outdoor walker, however is also compact enough for use indoors.

Features include a padded seat and padded folding backrest for comfort while seated. The adjustable height handles are fitted with Orthopaedic hand grips therefore providing the user with excellent support. The walker also folds for storage and transport. Two seat heights available.

## PRODUCT BENEFITS:

- Seat height of 55cm
- Height adjustable ergonomic shaped hand grips
- Lightweight aluminium construction
- Lockable handbrakes
- 200mm (8") castors

## SPECIFICATIONS:

TECHNICAL DATA	
Product Code	6865
Seat Width	440mm
Seat Height	550mm
Handle Height	790-910mm
Walker Weight	7.5kg
SWL	120kg

**activeMEDICAL**

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au  
activemedicalsupplies.com.au