

Clamp on Support Rail



For those needing support when standing from the bed, the clamp on rail is the solution. Rungs at correct spacing to prevent entrapment.

- 30cm wide x 40cm high - 10cm spacing between rungs
- Simply clamps to metal side frame of bed
- 150kg SWL
- Has been tested and conforms to AS/NS S3200.1.0 & AS/NZS 3200.2.38 Section 4
- Alternative size also available 7009 (30cm wide x 25cm high)

activeMEDICAL

SUPPLIES > SOLUTIONS > SUPPORT

35 Southgate Avenue Cannon Hill QLD 4170 | P +61 7 3899 8887 | F +61 7 3899 8825

enquiries@activemedicalsupplies.com.au
activemedicalsupplies.com.au